



THIS TEST BOX PLUGS INTO J2 ON THE REAR OF THE POWER MONITOR BOX.
 THE LED TURNS ON WHEN CURRENT IS FLOWING TO THE INDICATED LAMP

REVISIONS http://loel.ucolick.org/manual/Apf/schematics/APF_calib_lamp_sh_5.pdf
 02-23-11 REMOVED *2nd* HALOGEN LAMPS. NEW INTEGRATING SPHERE DESIGN REQUIRES A SINGLE (PLUS A LIVE SPARE) LAMP.

UNIVERSITY OF CALIFORNIA LICK OBSERVATORY		EL-3514 TEST BOX CALIBRATION LAMP STAGE APF SPECTROMETER	
DES'N BY: B. Alcott	ORIGIN DATE: 02-09-04	DWG. NO.	NUM. 5 OF 5
Last SN	MODIFY DATE: 02-23-11	EL-3514	REV. C
PRINT TIME: 06:37:45	PRINT DATE: Wed Feb 23, 2011		
PATH: APF\		APF_calib_lamp.sch	